





## Intraoral Auxiliaries

### INITIALLOY Springs

Made from the same alloy as INITIALLOY archwires, these springs exert constant and calibrated forces.



Load/flexion diagram for open springs



Load/flexion diagram for closed springs



#### INITIALLOY Stop Wound™ Open Springs

Each length has tight spirals every 4 cm. Once determining the measurement necessary, simply cut the wire at the center of the spiral. The spring goes back to its original size once cut. Springs are available in forces of 100 grams and 150 grams. Packaged by three, 18 cm length per tube.

Description	Force	Color	Internal bore	P/N
Extra-Light	50 g	■ Black	.036 "	60-0112-0000
Light	100 g	■ Blue	.036 "	60-0113-0000
Medium	150 g	■ Yellow	.036 "	60-0114-0000
Medium	150 g	■ Yellow	.045 "	60-0114-0045
Heavy	200 g	■ Red	.036 "	60-0115-0000



#### INITIALLOY Open Springs

These springs are 15 mm long before activated and can be compressed down to 3 mm, exerting a calibrated force until they recover their initial length. Springs are available in forces of 100 grams and 150 grams. Packaged by ten 15 cm lengths per tube.

Description	Force	Color	P/N
Light	100 g	■ Blue	60-0123-0000
Medium	150 g	■ Yellow	60-0124-0000